

## How to Use NasoNeb® Nasal Nebulizer




### **Step 1**

Arrange materials from package. **Carefully read instructions** on prescribed medication label before continuing. *Disclaimer: Your package may not contain some products appearing in this picture.*



**Step 2** Place the medication in the NasoNeb® cup as directed on prescription bottle.

**Step 3** Placing the nozzles in your nostrils, turn on machine and inhale medication. The NasoNeb® typically takes 1-2 minutes to deliver your prescribed medication.



### **Technique:**

Tilt your head slightly downward. Breathe through both your nose and mouth with the compressor on. Turn the compressor off when the main canister is empty.

### **Cleaning:**

Clean the nasal nebulizer with warm, soapy water after each use. Avoid the use of harsh detergents or cleaning agents. *Ensure the compressor is turned off and unplugged and product is disassembled.*

For more information on the NasoNeb® product, please visit: [www.mcpnc.com](http://www.mcpnc.com).

Medical Center Pharmacy  
912 S. 16<sup>th</sup> Street  
Wilmington, NC 28401  
(910)763-1896 / (800) 930-8785

