

## How to Use Nasal Atomizer



**Step 1:** Arrange materials from package. **Carefully read instructions** on prescribed medication label before continuing

**Step 2:** Add 10mL of medication to medicine cup and fill the syringe.




**Step 3:** Attach the nasal atomizer tip to the syringe by screwing on.

**Step 4.** Place the tip of the atomizer in the nostril. Press the plunger and gently spray half (5mL) of the medication into the nostril. Move the device over to the opposite nostril and repeat.

**Cleaning:**

1. Wash with warm, soapy water.
2. Rinse with clean water.
3. Allow to dry between each use.



*If you have any questions, please call Medical Center Pharmacy at 800.930.8785*

Medical Center Pharmacy  
912 S. 16<sup>th</sup> Street  
Wilmington, NC 28401  
(910)763-1896